



Putting your patient's life first

Aligned with the same platform as our System One sleep therapy technology, the Philips Respironics **BiPAP AVAPS noninvasive ventilator** provides the comfort of pressure ventilation and the consistent efficiency of assured tidal volume. Adapting automatically to changing patient needs, BiPAP AVAPS incorporates System One advanced technologies that help simplify patient care.

System One features

- Advanced event detection
- System One Resistance Control
- Oximetry module
- Encore and DirectView patient management systems
- SD card and wired, wireless, and broadband modem capable
- Digital Auto-Trak

AVAPS (Average Volume Assured Pressure Support)

- Automatically adapts to disease progressions and changing patient needs
- Maintains optimal patient comfort while supporting patient care and efficacy
- Improves ventilation efficacy
- Simplifies the titration process

Heated Tube humidification

- Air temperature control
- Improved humidification
- Rainout protection

System One

BiPAP AVAPS



Enhanced functionality

- Increased pressure capabilities to 30 cm H₂O
- Onboard memory stores patient therapy data without an SD card
- System One alignment leverages accessories across all Philips Respironics sleep therapy and noninvasive ventilator devices
- Heated Tube technology tailors humidity levels to patient needs

Home Healthcare
Respiratory Solutions

★ Company of the Year ★
2012-2014

awarded by
FROST & SULLIVAN


PHILIPS
RESPIRONICS

BiPAP AVAPS ordering information

Specifications		Device ordering information	Part number
Pressure range	EPAP 4-30 cm H ₂ O IPAP 4-30 cm H ₂ O CPAP 4-20 cm H ₂ O	BiPAP AVAPS, C Series	IN1161X
Breath rate	0-30 BPM	Ventilation accessories	
Flex pressure relief	0 to 3 (S-Mode)	Humidification	
Ramp time	0 to 45 min (5-min increments)	System One 60 Series humidifier	IN6H
Starting ramp	4 to EPAPmin	System One 60 Series Heated Tube humidifier (includes 80 watt power supply and Heated Tube)	IN6T
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H	Water chamber (includes both 50 and 60 Series lids)	1063785
Weight	3 lbs; 1.36 kg (without humidifier)	Heated Tube humidifier upgrade kit	6T-CONVKIT
Filters	Foam and optional ultra-fine	Tubing	
Device setup	LCD/control wheel/push button	System One Heated Tube	SYS1HT15
Data storage capacity (minimum)	Display: compliance VIC, 1-, 7- and 30-day averages SD card > 1 year On-board: up to 3 months	System One Performance Tubing (15mm)	1074400
Compliance meter	Breathing detection	System One Performance Tubing (22mm)	1032907
Altitude compensation	Automatic	Filters	
Electrical requirements	100-240 VAC, 50/60 Hz	Pollen (foam) filter, reusable (1/pack)	1035443
Humidification	Integrates with System One 60 Series humidifier and System One 60 Series Heated Tube humidifier	Pollen (foam) filter, reusable (2/pack)	1029330
Miscellaneous	Lighted LEDs	Ultra-fine filter, disposable (1/pack)	1063096
		Ultra-fine filter, disposable (2/pack)	1035442
		Ultra-fine filter, disposable (6/pack)	1029331
		Accessories	
		Device side panel	1099592
		Carrying case	1063857

*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).

Respironics, BiPAP, AVAPS, Bi-Flex, Encore and System One are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.



For more information, please visit www.philips.com/respironics



© 2013 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy AE 7/26/13 MCI 4105740 PN 1111568

Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics Benelux
+33 1 47 28 30 83
Philips Respironics Deutschland
+49 8152 93 06 0
Philips Respironics France
+33 2 51 89 36 00
Philips Respironics Iberia
+34 91 566 90 00

Philips Respironics Italy
+39 039 203 1
Philips Respironics Sweden
+46 8 120 45 900
Philips Respironics Switzerland
+41 6 27 45 17 50
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668 USA